

Pressure kit components & Use

Usage of the pressure kit

Pressure kit is enable the user to expend it use of the pressure sensor in a way that he/he can use any of the kit components to connect to any other device/s.

Pressure kit content in picture



Pressure kit components explanation

1. 13381 - 60 ml syringe with Luer Lock. The syringe can be used in order to perform experiments and also to calibrate the sensor.
2. 13882 - 50 cm transparent tubing 6 mm outer diameter.
3. 13883 - Syringe extender (x3).
4. 13880 - Three way valve (x2).
5. 14415 - male luer lock to barb for 1/8" tube.
6. 13879 - Stopper #5
7. 13873 - Stopper #3

